

PROSUM® Angiographic Syringe

Features & Benefits

- Crystal clear polycarbonate barrel with bold graduation provides clear and precise reading
- Smooth injections and aspirations
- Available with Rotating Adaptor provides easy and secure catheter connection
- Available with thumb ring or palm pad
- Available with 0.5 mL Reservoir - Prevents stopper from bottoming out during injection
- Supplied in sterile individual blister pack OR bulk non-sterile
- Other variants available upon request



PROSUM® Angiographic Syringe				
REF	Description			
	Size	Plunger Type	Fitting	0.5ml Reservoir
ANS060R1	6ml	Thumb Ring	Rotating adaptor	Yes
ANS060R0	6ml	Thumb Ring	Rotating adaptor	No
ANS060F1	6ml	Thumb Ring	Fixed MLL	Yes
ANS060F0	6ml	Thumb Ring	Fixed MLL	No
ANS080R1	8ml	Thumb Ring	Rotating adaptor	Yes
ANS080R0	8ml	Thumb Ring	Rotating adaptor	No
ANS080F1	8ml	Thumb Ring	Fixed MLL	Yes
ANS080F0	8ml	Thumb Ring	Fixed MLL	No
ANS100R1	10ml	Thumb Ring	Rotating adaptor	Yes
ANS100F1	10ml	Thumb Ring	Fixed MLL	Yes
ANS101R1	10ml	Palm Pad	Rotating adaptor	Yes
ANS101F1	10ml	Palm Pad	Fixed MLL	Yes
ANS120R1	12ml	Thumb Ring	Rotating adaptor	Yes
ANS120R0	12ml	Thumb Ring	Rotating adaptor	No
ANS120F1	12ml	Thumb Ring	Fixed MLL	Yes
ANS120F0	12ml	Thumb Ring	Fixed MLL	No
ANS121R1	12ml	Palm Pad	Rotating adaptor	Yes
ANS121F1	12ml	Palm Pad	Fixed MLL	Yes

Individual Blister Pack; 25/Box; 100/Carton

