

PROSUM® Inflation Device

Features & Benefits

- PROSUM® inflation device is a versatile device for a variety of interventional procedures such as balloon angioplasty
- Polycarbonate barrel
- ABS plunger
- Latex free
- Other variants available upon request



PROSUM® Inflation Device				
REF	Description			
	Volume	Max. Pressure	Tubing	Quick Release Button
IND022030	20ml	30 atm/bar	30cm High Pressure Tube with Male Rotating Luer Connector	YES
IND012030	20ml	30 atm/bar	30cm High Pressure Tube with Male Rotating Luer Connector	NO
IND012040	20ml	40 atm/bar	30cm High Pressure Tube with Male Rotating Luer Connector	NO





