







PROSUM[®] Insulin Syringe




Features & Benefits

- PROSUM[®] Insulin Syringes are available in a wide range of sizes and other variants
- Thin-walled needle with tricut bevel achieves almost pain free injection
- Transparent barrel for easy identification of air bubbles
- Bold graduation provides clear and precise reading
- Low dead space design for delivery accuracy and reduced medication waste
- Available in U-100 and U-40
- Latex free
- Supplied in sterile individual blister pack, tray pack, PE bag or bulk non-sterile
- Other variants available upon request

PROSUM [®] Insulin Syringe			
REF	Description		
0.3ml Insulin Syringe with Fixed Needle			
 U-100	 U-40		
INS31271271	INS34271271	0.3ml	27G × 1/2" (12.7mm)
INS31281271	INS34281271	0.3ml	28G × 1/2" (12.7mm)
INS31291271	INS34291271	0.3ml	29G × 1/2" (12.7mm)
INS31301271	INS34301271	0.3ml	30G × 1/2" (12.7mm)
INS31300801	INS34300801	0.3ml	30G × 5/16" (8mm)
INS31310801	INS34310801	0.3ml	31G × 5/16" (8mm)
0.5ml Insulin Syringe with Fixed Needle			
 U-100	 U-40		
INS51271271	INS54271271	0.5ml	27G × 1/2" (12.7mm)
INS51281271	INS54281271	0.5ml	28G × 1/2" (12.7mm)
INS51291271	INS54291271	0.5ml	29G × 1/2" (12.7mm)
INS51301271	INS54301271	0.5ml	30G × 1/2" (12.7mm)
INS51300801	INS54300801	0.5ml	30G × 5/16" (8mm)
INS51310801	INS54310801	0.5ml	31G × 5/16" (8mm)
1ml Insulin Syringe with Fixed Needle			
 U-100	 U-40		
INS11271271	INS14271271	1ml	27G × 1/2" (12.7mm)
INS11281271	INS14281271	1ml	28G × 1/2" (12.7mm)
INS11291271	INS14291271	1ml	29G × 1/2" (12.7mm)
INS11301271	INS14301271	1ml	30G × 1/2" (12.7mm)
INS11300801	INS14300801	1ml	30G × 5/16" (8mm)
INS11310801	INS14310801	1ml	31G × 5/16" (8mm)

Available in:
 Individual Blister Pack: 100/Box; 3200/Carton
 OR Tray Pack: 25/Tray; 1000/Carton
 OR PE bag: 10/PE bag; 100/Box; 3200/Carton



PROSUM [®] Insulin Syringe			
REF	Description		
1ml U-100 Insulin Syringe with Detachable Needle			
S1-01126050I	1ml		26G × 1/2" (12.7mm)
S1-01126580I	1ml		26G × 5/8" (16mm)
S1-01125580I	1ml		25G × 5/8" (16mm)
Individual Blister Pack: 100/Box; 3200/Carton			
1ml U-100 Insulin Syringe without Needle			
S0-010I	1ml		Luer lock
S0-011I	1ml		Luer Slip
Individual Blister Pack: 100/Box; 3600/Carton			



CE 0088 ~~LATEX~~ ~~PVC~~ ~~☒~~ STERILE EO