

PROSUM® IV Catheter

Features & Benefits

- PROSUM® IV Catheters are available in Straight, Winged and Ported variants for different applications
- Ultrasharp needle with precisely tapered catheter tip and the special lubrication ensure smooth and comfortable catheter advancement
- Catheters are available in polyurethane (PUR) for softer and more comfortable performance, or FEP with firmer construction for arterial access and special procedures
- Clear flashback Chamber allows immediate confirmation of venipuncture
- Removable flashplug allows syringe connection for aspiration or other special procedures
- Quick and easy gauge identification with colour coded hub
- DEHP-free, PVC-free and latex-free
- Other variants available upon request



Ultrasharp needle with special lubrication ensure smooth veinpuncture

Wing design for stabilization and easy handling

Integrated capped injection port for safe and easy medication injection without an extra needle or risk of needlestick injuries



PROSUM® Straight IV Catheter

REF		Description			Flow Rate
PUR	FEP	Size			
IVC114200P	IVC114200F	14G × 2"	(2.0 × 51mm)		310
IVC116200P	IVC116200F	16G × 2"	(1.7 × 51mm)		195
IVC118114P	IVC118114F	18G × 1 1/4"	(1.3 × 32mm)		103
IVC118112P	IVC118112F	18G × 1 1/2"	(1.3 × 38mm)		98
IVC118200P	IVC118200F	18G × 2"	(1.3 × 51mm)		85
IVC120114P	IVC120114F	20G × 1 1/4"	(1.1 × 32mm)		60
IVC120200P	IVC120200F	20G × 2"	(1.1 × 51mm)		53
IVC122100P	IVC122100F	22G × 1"	(0.9 × 25mm)		38
IVC124034P	IVC124034F	24G × 3/4"	(0.7 × 19mm)		22
IVC126035P	IVC126035F	26G × 3/5"	(0.6 × 15mm)		15

Individual Blister Pack; 50/Box; 500/Carton

PROSUM® Winged IV Catheter

REF		Description			Flow Rate
PUR	FEP	Size			
IVC214200P	IVC214200F	14G × 2"	(2.0 × 51mm)		310
IVC216200P	IVC216200F	16G × 2"	(1.7 × 51mm)		195
IVC218114P	IVC218114F	18G × 1 1/4"	(1.3 × 32mm)		103
IVC218112P	IVC218112F	18G × 1 1/2"	(1.3 × 38mm)		98
IVC218200P	IVC218200F	18G × 2"	(1.3 × 51mm)		85
IVC220114P	IVC220114F	20G × 1 1/4"	(1.1 × 32mm)		60
IVC220200P	IVC220200F	20G × 2"	(1.1 × 51mm)		53
IVC222100P	IVC222100F	22G × 1"	(0.9 × 25mm)		38
IVC224034P	IVC224034F	24G × 3/4"	(0.7 × 19mm)		22
IVC226035P	IVC226035F	26G × 3/5"	(0.6 × 15mm)		15

Individual Blister Pack; 50/Box; 500/Carton

PROSUM® Ported IV Catheter

REF		Description			Flow Rate
PUR	FEP	Size			
IVC314200P	IVC314200F	14G × 2"	(2.0 × 51mm)		310
IVC316200P	IVC316200F	16G × 2"	(1.7 × 51mm)		195
IVC318114P	IVC318114F	18G × 1 1/4"	(1.3 × 32mm)		103
IVC318112P	IVC318112F	18G × 1 1/2"	(1.3 × 38mm)		98
IVC318200P	IVC318200F	18G × 2"	(1.3 × 51mm)		85
IVC320114P	IVC320114F	20G × 1 1/4"	(1.1 × 32mm)		60
IVC320200P	IVC320200F	20G × 2"	(1.1 × 51mm)		53
IVC322100P	IVC322100F	22G × 1"	(0.9 × 25mm)		38
IVC324034P	IVC324034F	24G × 3/4"	(0.7 × 19mm)		22
IVC326035P	IVC326035F	26G × 3/5"	(0.6 × 15mm)		15

Individual Blister Pack; 50/Box; 500/Carton