

PROSUM®

Blunt Fill Needle & Blunt Fill Needle with Filter

Features & Benefits

- Ideal for use with glass vials/ampoules
- Significantly reduces risk of needlestick injury
- 10x the force to penetrate the skin, compared to a conventional needle
- 45 degree angle cut on bevel reduces coring of vial stopper
- 5 micron filter to capture particles during drawing up
- Bold pink coloured cap for easy identification from standard needles
- Sterile, Individual blister pack
- Other variants available upon request

PROSUM® Blunt Fill Needle	
REF	Description
BN20100	20G × 1" (0.9mm × 25mm)
BN19100	19G × 1" (1.1mm × 25mm)
BN18100	18G × 1" (1.2mm × 25mm)
BN18112	18G × 1 ½" (1.2mm × 40mm)
BN18200	18G × 2" (1.2mm × 50mm)

100/Box; 1000/Case



PROSUM® Blunt Fill Needle with Filter	
REF	Description
BNF20100	20G × 1" (0.9mm × 25mm) with 5 micron filter
BNF19100	19G × 1" (1.1mm × 25mm) with 5 micron filter
BNF18100	18G × 1" (1.2mm × 25mm) with 5 micron filter
BNF18112	18G × 1 ½" (1.2mm × 40mm) with 5 micron filter
BNF18200	18G × 2" (1.2mm × 50mm) with 5 micron filter

100/Box; 1000/Case






