

PROSUM® IV Accessories

Features & Benefits

- PROSUM® offers wide range of IV administration accessories to meet different medical practice
- All DEHP-free, Latex-free
- Supplied in sterile individual blister pack OR bulk non-sterile
- Other variants available upon request



IVA001
Dual Red Cap
Male & Female End
PP



IVA004Y
Injection Site (Heparin) Cap
Male Luer Lock & Injection Site
Yellow
PP, Polyisoprene



IVA004W
Injection Site (Heparin) Cap
Male Luer Lock & Injection Site
White
PC, Polyisoprene



IVA002B
Male Luer Cap, Blue, PP

IVA002W
Male Luer Cap, White, PP

IVA002N
Male Luer Cap, Natural, PP



IVA005
One-Way Valve
Female Luer to Male Luer
DEHP-free, Latex-free
PC, Silicone



IVA003B
Female Luer Cap, Blue, PP

IVA003W
Female Luer Cap, White, PP

IVA003N
Female Luer Cap, Natural, PP



EXT053002
Extension Tube with NeedleFree Y Connector
NeedleFree Y Connector, Injection Site Cap, Extension Line, Male Luer with Spin Lock, Slide Clamp, Female Luer Cap, 150mm



EXT053001
Extension Tube with NeedleFree Connector
NeedleFree Connector, Injection Site Cap, Extension Line, Male Luer with Spin Lock, Slide Clamp, Female Luer Cap, 150mm