

PROSUM® Infusion Set

Features & Benefits






- PROSUM® provides wide range of IV administration set to meet different clinical applications
- NeedleFree connector – full needlestick prevention and easy medicine administration
- Latex-free, DEHP-free, non-toxic
- Other variants available upon request



PROSUM® Infusion Set	
REF	Description
INF001001	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, Male luer lock, Tube: 3mm/4mm (ID/OD), 1400mm</p>
INF001002	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, Y injection site, Male luer lock, Tube: 3mm/4mm (ID/OD), 1500mm</p>
INF001003	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Check Valves, Roller clamp, Male luer lock, Tube: 3mm/4mm (ID/OD), 1500mm</p>
INF001004	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, Needle free Y injection site, Male luer lock, Tube: 3mm/4mm (ID/OD), 1500mm</p>
INF001005	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, 2 Needle free Y injection sites, Male luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INF001006	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, Check valve, 2 Needle free Y injection sites, Male luer lock, Tube: 3mm/4mm (ID/OD), 2000mm</p>
INF001007	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, 0.2 micron filter, Male luer lock, Tube: 3mm/4mm (ID/OD), 1400mm</p>
INF001008	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, Needle free Y injection site, 0.2 micron filter, Male luer lock, Tube: 3mm/4mm (ID/OD), 1500mm</p>
INF001009	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Roller clamp, 2 Needle free Y injection sites, 0.2 micron filter, Male luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>



50/Box; 400/Carton

PROSUM® Infusion Set (Continue)

REF	Description
INF002001	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Flow regulator, Slide clamp, Male spin luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INF002002	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Flow regulator, Slide clamp, Needle free Y injection site, Male spin luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INF002003	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Flow regulator, Slide clamp, 0.2 micron filter, Male spin luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INF002004	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Flow regulator, Slide clamp, 0.2 micron filter, Needle free Y injection site, Male spin luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INF002005	 <p>20 Drops/ml, Air filter, Drip chamber, Fluid filter, Flow regulator, Slide clamp, Check valve, 2 Needle free Y injection sites, Male spin luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>

50/Box; 400/Carton

PROSUM® Infusion Set with Burette

REF	Description
INFB10200	 <p>100ml burette, 20 Drops/ml, luer lock, Tube: 3mm/4mm (ID/OD), 1800mm</p>
INFB15200	150ml burette, 20 Drops/ml, luer lock, Tube: 3mm/4mm (ID/OD), 1800mm
INFB10600	100ml burette, 60 Drops/ml, luer lock, Tube: 3mm/4mm (ID/OD), 1800mm
INFB15600	150ml burette, 60 Drops/ml, luer lock, Tube: 3mm/4mm (ID/OD), 1800mm
INFB10201	 <p>100ml burette, 20 Drops/ml, Needle free Y injection site, luer lock, Tube: 3mm/4mm (ID/OD), 2000mm</p>
INFB15201	150ml burette, 20 Drops/ml, Needle free Y injection site, luer lock, Tube: 3mm/4mm (ID/OD), 2000mm
INFB10601	100ml burette, 60 Drops/ml, Needle free Y injection site, luer lock, Tube: 3mm/4mm (ID/OD), 2000mm
INFB15601	150ml burette, 60 Drops/ml, Needle free Y injection site, luer lock, Tube: 3mm/4mm (ID/OD), 2000mm

25/Box; 200/Carton