

PROSUM® Extension Set

Features & Benefits


- PROSUM® offers wide range of extension sets and accessories to meet different clinical applications
- NeedleFree connector – full needlestick prevention and easy medicine administration
- Latex-free, DEHP-free, non-toxic
- Other variants available upon request

CE 0088 ~~LATEX~~  

PROSUM® Extension Set			
REF	Description		
	Tube ID/OD	Tube Length	Feature
EXT011001	1.3mm/2.4mm	450mm	Female Luer, Male Luer
EXT011002		600mm	
EXT011003		800mm	
EXT011005		1500mm	
EXT011007		2000mm	
EXT012001	1.5mm/2.7mm	450mm	
EXT012002		600mm	
EXT012003		800mm	
EXT012005		1500mm	
EXT012007		2000mm	
EXT013001	3mm/4mm	250mm	
EXT013002		500mm	
EXT013003		1000mm	
EXT013004		1500mm	
100/Box; 400/Carton			
EXT021001	3mm/4mm	100mm	Stopcock, Male Luer
EXT021002		250mm	
EXT021003		300mm	
EXT021004		500mm	
EXT021005		800mm	
EXT021006		1000mm	
EXT021008		1500mm	
50/Box; 200/Carton			
EXT032001	3mm/4mm	150mm	Female Luer, Needle free 'Y' connector, Male Luer
EXT032002	3mm/4mm	250mm	Female Luer, 2 × Needle free 'Y' connector, Male Luer
100/Box; 400/Carton			


PROSUM® Extension Set (Continue)

REF	Description		
	Tube ID/OD	Tube Length	Feature
EXT052004	1.3mm/ 2.4mm	100mm	Needle free connector, Male Luer
EXT052001		100mm	2 × Needle free connector, Male Luer
EXT052002		100mm	3 × Needle free connector, Male Luer
EXT052003		100mm	4 × Needle free connector, Male Luer
EXT051001	3mm/4mm	100mm	Needle free connector, Male Luer
EXT051002		100mm	2 × Needle free connector, Male Luer
EXT051003		100mm	3 × Needle free connector, Male Luer
EXT051004		100mm	4 × Needle free connector, Male Luer




100/Box; 400/Carton

EXT061001	3mm/4mm	100mm	2 Port Manifold, Male Luer
EXT062001			3 Port Manifold, Male Luer
EXT063001			4 Port Manifold, Male Luer




50/Box; 200/Carton

EXT081001	3mm/4mm	100mm	Female Luer, Check Valve, Male Luer
EXT091001	3mm/4mm	100mm	Female Luer, 0.2 Micron Filter, Male Luer



100/Box; 400/Carton

EXT101001	3mm/4mm	100mm	Female Luer, Flow Regulator, Needle free Y-site, Male Luer
EXT101002	3mm/4mm	100mm	Female Luer, Flow Regulator, Y injection site, Male Luer



50/Box; 200/Carton

