


# PROSUM® Dental Needle

### Features & Benefits

- Excellent tissue penetration achieved with a high quality triple-bevel needle with polished and siliconized cannula
- High-grade surgical stainless steel tubing to reduce risk of breakage
- Bevel mark to ensure good orientation of the bevel
- Colour coded label for easy size identification
- Other variants available upon request



PROSUM® Dental Needle		
REF	Description	
DN03012	30G × 12mm	Extra Short
DN03020	30G × 20mm	Short
DN03025	30G × 25mm	Long
DN02725	27G × 25mm	Short
DN02732	27G × 32mm	Long
DN02738	27G × 38mm	Extra Long
DN02525	25G × 25mm	Short
DN02532	25G × 32mm	Long
DN02538	25G × 38mm	Extra Long

*100/Box; 5000/Carton*