

PROSUM® Syringe with Needle - 0.5ml & 1ml

Features & Benefits

- PROSUM® Syringes are available in a wide range of sizes and variants
- Ultra sharp needle minimizes injection pain and increase patient comfort
- Colour coded hubs for easy gauge identification
- Transparent barrel and bold graduation provides clear and precise reading
- Special lubrication offers smooth plunger movement
- Low dead space design for delivery accuracy and reduced medication waste
- Latex free
- Supplied in sterile individual blister pack OR bulk non-sterile
- Other variants available upon request



Low Dead Space Design

PROSUM® Tuberculin Syringe

REF	Description	
0.5ml Tuberculin Syringe with Fixed Needle		
SF-00527050T	0.5ml	27G × 1/2"
SF-00527580T	0.5ml	27G × 5/8"
1ml Tuberculin Syringe (Luer Slip) with Detachable Needle		
S1-01127050T	1ml	27G × 1/2"
S1-01127580T	1ml	27G × 5/8"
S1-01126050T	1ml	26G × 1/2"
S1-01126580T	1ml	26G × 5/8"
S1-01125580T	1ml	25G × 5/8"
1ml Tuberculin Syringe without Needle		
S0-010T	1ml	Luer lock
S0-011T	1ml	Luer Slip

Individual Blister Pack: 100/Box; 3200/Carton
Tray Pack: 25/Tray; 1000/Carton



PROSUM® Allergy Syringe

REF	Description	
0.5ml Allergy Syringe with Fixed Needle		
SF-00527050A	0.5ml	27G × 1/2"
SF-00527580A	0.5ml	27G × 5/8"
1ml Allergy Syringe with Fixed Needle		
SF-01028050A	1ml	27G × 1/2"
SF-01027050A	1ml	27G × 5/8"
SF-01026050A	1ml	26G × 1/2"
SF-01026580A	1ml	26G × 5/8"

Individual Blister Pack: 100/Box; 3200/Carton
Tray Pack: 25/Tray; 1000/Carton

PROSUM® 1ml Syringe with Detachable Needle

REF		Description	
Luer Lock	Luer Slip	Size	Needle
S1-01027050	S1-01127050	1ml	27G × 1/2"
S1-01027580	S1-01127580	1ml	27G × 5/8"
S1-01026050	S1-01126050	1ml	26G × 1/2"
S1-01026580	S1-01126580	1ml	26G × 5/8"
S1-01025580	S1-01125580	1ml	25G × 5/8"

Individual Blister Pack; 100/Box; 3200/Carton

