

# PROSUM® Polycarbonate & Vacuum Lock Syringe

### Features & Benefits

- PROSUM® Polycarbonate Syringe is ideal for special clinical applications such as critical care, cosmetic application, etc.
- Crystal clear polycarbonate barrel improves visibility of fluid level
- Bold graduation provides clear and precise reading
- Luer lock connector assure secure connection to needles or IV ports
- ABS plunger, available in various colours
- Latex-free, PVC-free, non-toxic
- Supplied in sterile individual blister pack OR bulk non-sterile
- Other variants available upon request



### PROSUM® Polycarbonate Syringe

REF	Description	Qty
SPC0-010	1ml Luer lock	100/Box; 3600/Carton
SPC0-030	3ml Luer lock	100/Box; 3200/Carton
SPC0-050	5ml Luer lock	100/Box; 2400/Carton
SPC0-100	10ml Luer lock	100/Box; 1500/Carton
SPC0-200	20ml Luer lock	100/Box; 900/Carton
SPC0-300	30ml Luer lock	100/Box; 600/Carton

Individual Blister Pack


### PROSUM® Vacuum Lock Syringe

REF	Description	Qty
SPCL0-100	10ml Luer lock	100/Box; 1500/Carton
SPCL0-200	20ml Luer lock	100/Box; 900/Carton
SPCL0-300	30ml Luer lock	100/Box; 600/Carton

Individual Blister Pack

### Vacuum Lock Syringe

built with locking mechanism to maintain negative pressure.



**Fixed Male Luer Lock**  
ensure secure connection

**Silicone Plunger Tip**

**Polycarbonate Barrel**  
for clear visibility

**ABS Plunger**  
with locking positions for fixing vacuum at different volume levels

